



Combined antimony electrode

6.0421.100

Combined Sb electrode. This electrode is well suited for pH titrations in:

- matrices containing hydrofluoric acid
- strongly hygroscopic matrices

Thanks to the robust/break-proof plastic shaft made of polypropylene, this electrode is mechanically resistant.

Phạm vi giao hàng 6.0421.100

Qt.	Order no.	Description
-----	-----------	-------------

100 Stk 6.1236.050 Sleeve with SGJ 14/12 mm
Sleeve with SGJ 14/12 mm, polyethylene.



100 Stk 6.2008.040 Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



100 Stk 6.2802.000 Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 µm)



Phụ kiện tùy chọn

Order no.	Description	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	