




Equipment for pH Adjustment at the 845 Eluent Synthesizer

6.5329.000

pH adjustment to Eluent Synthesizer. Includes pH electrode and additional accessories.

Phạm vi giao hàng 6.5329.000

Order Number	Description	Quantity
6.1446.000	SGJ stopper / B-14/(15)	1 PCS



6.1446.040	Threaded stopper / M6	2 PCS
------------	-----------------------	-------



6.1602.190	Eluent Bottle Head, 1 x SGJ 14, 4 x M6	1 PCS
------------	--	-------




Eluent bottle head to connect up to 4 tubings with M6 connectors. For eluent bottle of the Eluent synthesizer.

6.1838.010

Connecting Tubing to 845, 85 cm

1 PCS

Connection eluent bottle at the eluent Synthesizer - Titrino



6.2148.010

Remote Box MSB

1 PCS

Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.



6.2726.090

Antidiffusion valve

1 PCS

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140, 6.1543.150, 6.1543.160, 6.1543.200



6.3026.210

Exchange Unit 10 mL

1 PCS

Exchange unit with integrated data chip with 10 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

