



CVS electrode equipment with 2 mm platinum electrode for Professional VA/CVS instruments

6.5339.010

Complete electrode set for the determination of organic additives in electroplating baths with CVS (Cyclic Voltammetric Stripping). Contains rotating platinum working electrode (electrode diameter: 2 mm), reference electrode, platinum auxiliary electrode, measuring vessels and electrolyte solutions.

Phạm vi giao hàng 6.5339.010

Order Number	Description	Quantity
6.0343.100	Separate Pt Rod Electrode Platinum auxiliary electrode for voltammetry and CVS.	1 PCS
		
6.0728.130	LL-Ag/AgCl reference electrode for CVS Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.	1 PCS
		
6.1204.610	2 mm platinum electrode tip for CVS Metal electrode tip made of platinum (Pt), diameter of electrode disk $2 \text{ mm} \pm 0.05 \text{ mm}$, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.3 mm. Temperature range: 0 - 50 °C.	1 PCS
		

6.1204.510

Driving axle for rotating disk electrode (RDE)

1 PCS

Driving axle for rotating disk electrode (RDE) with titanium axle for Professional VA/CVS instruments, M3 thread.



6.1245.010

Electrolyte vessel made of plastic

1 PCS

Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.



6.1415.210

Titration vessel / 10-90 mL

1 PCS



6.1415.250

Titration vessel / 50-150 mL

1 PCS



6.2308.020

Electrolyte 3 mol/L KCl (250 mL)

1 PCS

Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)



6.2310.010

Electrolyte KNO₃-1mol/L 250 mL

1 PCS

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)

