



VA accessory equipment with SPE electrode shaft for professional VA instruments

6.5339.070

Accessory equipment for use with *screen-printed electrodes* (SPE). Contains electrode shaft for screen-printed electrodes, stirrer, and measuring vessel. Without electrodes.

Phạm vi giao hàng 6.5339.070

Order Number	Description	Quantity
6.1204.500	Stirrer for Professional VA instruments Stirrer for Professional VA instruments, complete, metal-free.	1 PCS
		
6.1244.060	Silicone seal for SPE electrode shaft (10 pieces) Seal made of silicone for SPE electrode shaft 6.1241.090 and measuring head (SPE) 6.01256.020, 10 pieces.	1 PCS
		
6.1241.090	Electrode shaft for screen-printed electrodes (SPE) Holder for <i>screen-printed electrodes</i> (SPE) with plug contacts for connecting the electrode cables for working, reference, and auxiliary electrodes. When used with screen-printed electrodes, the electrode shaft provides a complete sensor for voltammetry. Whereas the screen-printed electrodes must be replaced at regular intervals, the electrode shaft can continue to be used.	1 PCS
		

6.1415.150

Titration vessel / 5-70 mL

1 PCS

