





Accessories for 2.771.0210

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.771.0210

Qty.	Order no.	Description	
1 pcs	1.771.0010	771 IC Compact Interface	 <p>The universally applicable analog-digital converter with two channels and 24 Bit resolution. The data is passed along to the PC through the RS-232C interface and can be processed further there with IC Net or with MagIC Net.</p>
2 pcs	6.2128.170	Cable 771 - analog	 <p>For connecting 771 IC Compact Interface to analog output of the detector (732, 791)</p>
1 pcs	6.2134.100	RS-232 cable, 5 m, DB9 – DB9	 <p>Connecting Metrohm RS-232 interface to Titrinos, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.</p>
1 pcs	6.6059.321	MagIC Net 3.2 Compact CD: 1 license	 <p>Professional PC program for controlling one intelligent Compact IC instrument, one detector and one Autosampler, or one 771 Compact Interface. The software permits the inspection, data acquisition, data evaluation and data monitoring as well as report generation of ion chromatographic analyses.</p> <p>Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.</p> <p>MagIC Net Compact complies fully with FDA Regulation 21 CFR Part 11 as well as GLP.</p> <p>MagIC Net is available in 16 dialog languages.</p> <p>1 license.</p>

Optional accessories

Order no.	Description	
6.2134.140	Cable 771 (analog) - 819 Connection of analog and start signals to 771	
6.2134.160	Cable Lambda 1010 -762/830-819 Cable Lambda 1010 - IC Interface - IC Detector	
6.2134.190	Analog Cable Compact IC Instruments - Compact IC Interface Connecting cable Compact IC (Conductivity and UV/VIS) - Compact IC Interface	
6.2156.050	Analog cable IC Detector - Bioscan - IC Interface Connecting cable Bioscan - IC Interface for use with MagIC Net™ software	