





Accessories for 2.775.0130

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.775.0130

Qty.	Order no.	Description	
1 pcs	1.728.0010	728 Stirrer without stand For mounting on a stand rod 10 mm dia.; without base plate (or holder) and without stand rod; without power adapter, with fixed connecting cable.	
1 pcs	1.775.0010	775 Dosimat Digital buret for manual titrations and dispensing, with automatic stopcock switching, including push-button cable for manual dispensing control and 6.3026.220 Exchange Unit 20 mL.	
1 pcs	6.1903.030	Stirring bar / 25 mm Stirring bar with magnetic core, PTFE covering, length 25 mm.	
1 pcs	6.2001.010	Support stand with support rod For Dosimats and Titrinos	

1 pcs **6.2013.010** **Clamping ring**

For support rods with a diameter of 10 mm.



1 pcs **6.2021.020** **Electrode holder**

Electrode holder for 4 electrodes and 2 buret tips



1 pcs **6.2107.000** **Push-button cable for Dosimat**

For manual operation of Dosimats



1 pcs **6.2739.010** **Wrench**

For exchange units.



1 pcs **6.3026.220** **Exchange Unit 20 mL**

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE/PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

