





Accessories for 2.823.0030

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.823.0030

Qty.	Order no.	Description	
1 pcs	1.823.0010	823 Membrane Pump Unit Membrane pump for direct connection to Sample Processors.	
1 pcs	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	
2 pcs	6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
3 pcs	6.1805.060	FEP tubing / M6 / 60 cm With light and kink protection	
1 pcs	6.1808.170	Distributor for Sample Processors for rinsing and aspirating Distributor for installation on Sample Processors for rinsing and aspirating when using externally connected pumps.	

1 pcs **6.1812.000** **PTFE tubing 4 mm / 6 mm / 4 m**



1 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.



3 pcs **6.2740.020** **Spray nozzle**

For the fine-spraying of the rinsing solution.

