






Accessories for 2.827.0114

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.827.0114

| Qty. | Order no. | Description | |
|-------|------------|--|---|
| 1 pcs | 6.0228.010 | Primatrode with NTC (fixed cable 1.2 m) |  |
| | | <p>Combined pH electrode with integrated temperature sensor (NTC) and fixed cable (1.2 m). This electrode is suitable for an introduction to GLP-compliant pH measurements in solutions that do not contain precipitates, proteins, or sulfides. This electrode is mechanically resistant thanks to the robust/unbreakable plastic shaft made of polypropylene and impact protection for the glass membrane. Reference electrolyte: c(KCl) = 3 mol/L, storage in storage solution.</p> | |
| 1 pcs | 6.2013.010 | Clamping ring |  |
| | | For support rods with a diameter of 10 mm. | |
| 1 pcs | 6.2016.050 | Support rod / 300 mm |  |
| 1 pcs | 6.2021.020 | Electrode holder |  |
| | | Electrode holder for 4 electrodes and 2 buret tips | |
| 1 pcs | 6.2161.010 | Plug-in power supply unit EU 230V /6V DC |  |
| | | Power supply for 827 pH lab | |



1 pcs

6.2621.070

5 mm hex key for IC Sample Processors



Optional accessories

| Order no. | Description | |
|-------------------|---|---|
| 6.2307.230 | Buffer solutions pH 4, 7 and 9 Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL |  |
| 6.2325.000 | pHit kit Maintenance kit for pH electrodes The kit contains: <ul style="list-style-type: none">• 50 mL cleaning solution• 50 mL 3M KCl solution• 50 mL storage solution• 2 Storage vessels• Instructions for use |  |