


Accessories for 2.838.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.838.0020

Qty.	Order no.	Description	
1 pcs	1.838.0020	838 Advanced IC Injection Sample Processor Advanced IC Sample Processor with six-way injector for applications that require an additional injector.	
3 pcs	6.1608.080	Bottle / 300 mL Rinsing vessel for IC Sample Processors.	
1 pcs	6.1825.210	PEEK sample loop 20 µL For injection valve, with 2 PEEK pressure screws	
2 pcs	6.1826.360	Pump tubing LFL (white/white), 3 stoppers For Sample Processors	

1 pcs **6.1831.010** **PEEK capillary 0.25 mm i.d. / 3 m**

For all IC components.



1 pcs **6.1831.050** **PEEK capillary 0.5 mm i.d., 40 cm**

For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler



2 pcs **6.1831.060** **PEEK capillary 0.5 mm i.d., 1 m**

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



1 pcs **6.2064.000** **Beaker retainer for Advanced IC Sample Processor**

Retaining plate for 838 Advanced Sample Processors. For holding 2.5 or 11 mL vials.



1 pcs **6.2134.100** **RS-232 cable, 5 m, DB9 – DB9**

Connecting Metrohm RS-232 interface to Titrinos, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



1 pcs **6.2141.190** **Cable 732/819-766/788/838**
Cable IC Detector - Advanced Sample Processor



1 pcs **6.2621.300** **Tool set for Advanced IC Sample Processors**
Tool set consisting of open-end wrenches and hexagon keys.



1 pcs **6.2739.000** **Wrench**
For tightening connectors



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



1 pcs **6.2744.110** **Pump tubing - inline filter PEEK adapter**
Connector between pump tubing and inline filter



1 pcs **6.2744.160** **Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



1 pcs **6.2751.110** **Safety shield**

Safety shield for the Advanced IC/VA Sample Processors.






1 pcs **6.2846.010** **Sample needle made of zirconium oxide with PEEK tip**

Sample needle for Sample Processors and Autosampler in ion chromatography. The PEEK tip guarantees secure piercing of the vial covering.



Optional accessories

Order no.	Description	
6.1562.130	Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors	
	Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.	
6.2041.310	Sample rack 12 x 250 mL	
	Sample rack for 12 x 250 mL, one row	
6.2041.320	Sample rack 16 x 150 mL	
	Sample rack for 16 x 150 mL samples, one-row	
6.2041.340	Sample rack 24 x 75 mL	
	Sample rack for 24 x 75 mL samples, one-row.	
6.2041.350	Sample rack 48 x 75 mL	
	Sample rack for 48 x 75 mL samples, two-row	
6.2041.360	Sample rack for 12 x 150 mL	
	Sample rack for 12 x 150 mL samples, one-row	
6.2041.370	Sample rack 14 x 200 mL	
	Sample rack for 14 x 200 mL samples, one-row	

6.2041.380 Sample rack 14 x 235 mL

Sample rack for 14 x 235 mL samples, one-row



6.2041.430 Sample rack 127 x 11 mL

Sample rack for 127 x 11 mL samples with 2 rinsing beakers.



6.2041.440 Sample rack 148 x 11 mL + 3 x 300 mL

Sample rack for 148 x 11 mL samples with 3 rinsing beakers.



6.2041.450 Sample rack 56 x 11 mL + 56 x 50 mL

Multi-row sample rack for 56 samples in 11 mL sample vessels (6.2743.050) and 56 samples in sample vessels of 50 mL (6.2747.010).



6.2041.760 Sample rack 54 x 11 mL + 1 x 300 mL

Sample rack with small diameter for 54 x 11 mL sample vessels in two rows and additionally 1 x 300 mL bottle.



6.2629.000 Screw nut to needle holder IC

On Sample Processors



6.2743.040 Sample tubes 2.5 mL

2.5 mL. 2000 pieces. For IC Sample Processors



6.2743.047 Sample tubes 2.5 mL

2.5 mL. 200 pieces. For IC Sample Processors



6.2743.050 Sample tubes 11 mL

2000 pieces. For IC Sample Processors and VA Autosampler



6.2743.057 Sample tubes 11 mL

200 pieces. For IC Sample Processors and VA Autosampler



6.2743.070 Stopper with perforation

2000 pieces. For IC Sample Processors



6.2743.077 Stopper with perforation

200 pieces. For IC Sample Processors

