





Accessories for 2.844.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.844.0010

Qty.	Order no.	Description	
1 pcs	1.844.0010	844 UV/VIS Compact IC Compact ion chromatograph for the photometric determination of polar substances, anions and cations in the UV/VIS range. Very compact design. System with integrated dual-piston high-pressure pump, electrically operated injection valve and highly sensitive diode array detector.	
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
1 pcs	6.1608.070	Eluent bottle / 2 L / GL 45	
1 pcs	6.1609.000	Adsorbing tube, large and bent For filling with adsorber material.	

1 pcs **6.1816.060** **Drain tubing for UV/VIS Compact IC and Bioscan**
Drainage tubing for draining off escaped liquid



1 pcs **6.1831.010** **PEEK capillary 0.25 mm i.d. / 3 m**
For all IC components.



2 pcs **6.1831.110** **PEEK capillary 0.25 mm i.d., 30 cm**



1 pcs **6.1834.010** **PTFE aspiration tube to IC Pump, 2.5 m**
For Compact IC, Personal IC, Basic IC, 6.5324.000 bottle holder.



1 pcs **6.2023.020** **Clip for SGJ 14/15**
Clip for SGJ 14/16



1 pcs **6.2027.030** **Column holder 8.5 mm**
For mounting columns in column chamber of IC instruments.



1 pcs **6.2027.040** **Column holder 11.3 mm**
For mounting columns in column chamber of IC instruments.



1 pcs **6.2134.100** **RS-232 cable, 5 m, DB9 – DB9**
Connecting Metrohm RS-232 interface to Titrinos, 756 KF Coulometer, 762 or 771 IC Interfaces, 790 Personal IC, 743 Rancimat, 763 PVC Thermomat, 873 Biodiesel Rancimat.



1 pcs **6.2617.010** **Tool for piston seal**
For removing and assembling the piston seal for all standard pump heads



1 pcs **6.2620.150** **Pulsation absorber (metal-free)**
For all ion chromatographs



1 pcs **6.2621.000** **Adjustable wrench**
Max. opening: 20 mm. For IC instruments



1 pcs **6.2621.030** **Hexagon key 4 mm**



1 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



1 pcs **6.2816.020** **Syringe 10 mL with Luer connection**
For various applications in IC and VA



1 pcs **6.2821.090** **Aspiration filter**
Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.




1 pcs **6.2821.130** **Spare filter for inline filter**
Spare filter for inline filter.



1 pcs **6.2839.130** **Flow cell 10 mm to UV/VIS instruments**
PEEK flow cell



Optional accessories

Order no.	Description	
6.2804.060	Deuterium lamp (UV) for UV/VIS instruments Replacement bulb for the UV range	
6.2839.100	Flow-through cell (stainless steel) to UV/VIS Compact IC	