


## Accessories for 2.852.0160

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.852.0160

Qty.	Order no.	Description	
1 pcs	1.852.0010	852 Titrandos	
1 pcs	1.900.0010	900 Touch Control	
		Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick.	
		Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian.	
1 pcs	6.0338.100	Double Pt-wire electrode for volumetry	
		Indicator electrode used for volumetric Karl Fischer titration.	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry	
		Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0345.100	Generator electrode without diaphragm	
		Generator electrode without diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for most applications where the sample is highly soluble in alcohol.	

---

**1 pcs    6.1244.040    Set of seals**

Set of seals for 6.1414.030. The set is comprised of one seal and various O-rings.



---

**1 pcs    6.1403.030    KF adsorber tube for coulometer cell**

For KF coulometers.



---

**1 pcs    6.1403.040    KF Adsorber tube**

Used with KF equipment. Including cover and O-ring.



---




**1 pcs    6.1414.030    KF titration vessel lid**

For volumetric KF titrations.



---

**1 pcs    6.1415.220    Titration vessel / 20-90 mL**

1 pcs	6.1415.250	Titration vessel / 50-150 mL	
1 pcs	6.1437.000	SGJ stopper SGJ 14	
1 pcs	6.1446.060	Stopper / B-14/15 / M10	
2 pcs	6.1448.010 Set of 5 items.	Septum 12 mm, 5 pieces	
2 pcs	6.1448.020 Set of 5 pcs.	Septum 16 mm, 5 pieces	
1 pcs	6.1464.320	KF titration vessel / 80-250 mL / coulometric	

---

**2 pcs**     **6.1903.020**     **Stirring bar / 16 mm**  
Stirring Bar with magnetic core, PTFE coated, length 16 mm



---

**3 pcs**     **6.1903.030**     **Stirring bar / 25 mm**  
Stirring bar with magnetic core, PTFE covering, length 25 mm.



---

**1 pcs**     **6.2001.090**     **Base plate on both sides**  
Double stand for attaching 2 x 801, 803, 804 to a Titrande



---

**2 pcs**     **6.2013.010**     **Clamping ring**  
For support rods with a diameter of 10 mm.



---

**2 pcs**      **6.2016.030**      **Support rod / 250 mm**



---

**2 pcs**      **6.2043.005**      **Holding clip for bottles**

Holding spring for reagent bottles in exchange units.



---

**1 pcs**      **6.2047.020**      **Titration vessel holder**

Titration vessel holder for coulometric cells. Used with 756, 831 KF Coulometer



---

**2 pcs**      **6.2103.130**      **Adapter red 2 mm plug / 4 mm socket**

For connecting plug B (4 mm) to 2 mm socket.



---

**2 pcs**      **6.2103.140**      **Adapter black 2 mm plug / B socket 4 mm**

For connecting plug B (4 mm) to 2 mm socket.



2 pcs	6.2104.020	<b>Electrode cable / 1 m / F</b> For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
1 pcs	6.2104.120	<b>Electrode cable / 1 m / H</b> With plug H. For connecting generator electrode - KF Coulometer.	
1 pcs	6.2412.000	<b>Weighing boats</b>  Weighing boat made of glass for the 6.1414.030 titration vessel lid or the 6.1455.31X, 6.1464.32X, 6.1465.320 titration vessel. For KF water determinations. Including protective tube.	
1 pcs	6.2701.040	<b>GL 18 screw cap with hole</b> Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.	
3 pcs	6.2713.000	<b>Sleeve</b>  For the fat-free ground-joint connection with the standard ground-joint SGJ 14.	
1 pcs	6.2713.010	<b>Sleeve</b>  Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.	

---

**1 pcs**    **6.2713.020**    **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



---

**2 pcs**    **6.2730.010**    **Screw nipple**

Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.



---

**1 pcs**    **6.2730.020**    **Septum stopper**

With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.



---

**4 pcs**    **6.2730.030**    **Stopper**

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



---

**1 pcs**    **6.2738.000**    **Funnel**

For the KF coulometer.



---

**1 pcs    6.2811.000    Molecular sieve**

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



---

**2 pcs    6.2816.030    Needle with Luer connector**

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



---

**2 pcs    6.2816.090    Syringe 5 mL**

5 mL syringe made of PP, with Luer connector.



---

**1 pcs    6.6064.010    USB Memory Stick for 900 Touch Control**