






Accessories for 2.863.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.863.0020

Qty.	Order no.	Description	
1 pcs	1.863.0010	863 Compact IC Autosampler The 863 Compact IC Autosampler is the ideal help for routine analysis. It enables automation of the full range of Metrohm IC instruments.	
1 pcs	6.1822.410	PEEK connection capillary Connection capillary for sample transfer with sample changers, 1 x M6, inner diameter 0.5 mm, for trace analysis	
2 pcs	6.1826.020	PVC pump tubing (blue/blue) For online IC instruments and automation in voltammetry.	
1 pcs	6.1831.050	PEEK capillary 0.5 mm i.d., 40 cm For inline dialysis and for sample changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit, Compact VA Autosampler	
1 pcs	6.2141.230	Remote cable Compact Sample Changer/Compact Autosampler - Titrino plus/Pump Station Cable for the remote connection between Compact Sample Changer/Compact Autosampler and Titrino plus and/or Pump Station.	

1 pcs **6.2151.100 Adapter USB MINI (OTG) - USB A**
For connecting USB instruments.



1 pcs **6.2621.090 1/2 in. wrench**
For 1/2 in. screws. For IC Sample Processors



1 pcs **6.2621.130 Hexagon key 2 mm**
2 mm.



1 pcs **6.2621.140 Hexagon key 2.5 mm**



1 pcs **6.2739.000 Wrench**
For tightening connectors



1 pcs **6.2743.057 Sample tubes 11 mL**
200 pieces. For IC Sample Processors and VA Autosampler



1 pcs **6.2744.070 Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



2 pcs **6.2744.160 Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



1 pcs **6.2833.030 Needle holder 1/8 in. for 863 Compact Autosampler**
Needle holder 1/8 in. for 863 Compact Autosampler



1 pcs 6.2846.000 Sample needle made of zirconium oxide

Sample needle for Sample Processors and Autosampler in ion chromatography and voltammetry.



Optional accessories

Order no.	Description	
2.843.0040	843 VA Membrane Pump Station for Compact VA systems	
	<p>The 843 Membrane Pump Station has two built-in membrane pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace.</p> <p>The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.</p>	
2.843.0140	843 VA Peristaltic Pump Station for Compact VA systems	
	<p>The 843 Peristaltic Pump Station has two built-in peristaltic pumps that can be controlled directly by the 797 VA Computrace via the remote signal. It is used for the automatic rinsing and emptying of the measuring vessel in partly automated and fully automated VA and CVS systems with the 797 VA Computrace.</p> <p>The version for voltammetry is equipped with the complete accessories for the automatic rinsing and emptying of the measuring vessel.</p>	
6.2723.310	Protective cover for Titrino plus	
	<p>Protective cover for display and keypad for all instruments of the Titrino plus family.</p> <p>Suitable for: 848 Titrino plus, 862 Compact Titrosampler, 863 Compact Autosampler, 865 Dosimat plus, 869 Compact Sample Changer, 870 KF Titrino plus, 876 Dosimat plus, 877 Titrino plus, 899 Coulometer.</p>	
6.2743.050	Sample tubes 11 mL	
	<p>2000 pieces. For IC Sample Processors and VA Autosampler</p>	