

Accessories for 2.875.9020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.



Scope of delivery 2.875.9020

Qty.	Order no.	Description
1 pcs	1.875.9020	KF Gas Analyzer with TFT monitor

The 875 KF Gas Analyzer is a completely configured analysis system. The instrument contains a computer with **tiamo**™ full titration software installed, a monitor with integrated keyboard and a mouse pad, in addition to an 851 Titrando for automatic coulometric water content determination. Predefined methods for water content determination in liquefied and permanent gases facilitates work with the 875 KF Gas Analyzer. In addition to the determination of water content in gases, the system can also be used for the determination of water content in liquid and solid samples.



Optional accessories

Order no.	Description	
6.7209.000	KF Gas Analyzer equipment for reagent replacement and adding methanol	
	KF Gas Analyzer accessories for automatic reagent replacement or the addition of methanol.	
6.7209.010	KF Gas Analyzer equipment for rinsing with solvent	
	KF Gas Analyzer accessories for rinsing the lines with solvent.	