






Accessories for 2.900.0010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.900.0010

Qty.	Order no.	Description	
1 pcs	1.900.0010	<p>900 Touch Control</p> <p>Operating unit for the Titrandos, USB Sample Processors, 856 Conductivity Module, 867 pH Module and 846 Dosing Interface. Touch-sensitive, high-resolution color display, simple and intuitive operation, thanks to Favorites for direct method access. With integrated Ethernet interface for direct connection to the Internet and USB interface for connecting USB printers or a USB memory stick.</p> <p>Dialog languages: German, English, Chinese, French, Spanish, Portuguese, Russian, Korean, Polish and Italian.</p>	
1 pcs	6.6064.010	<p>USB Memory Stick for 900 Touch Control</p>	

Optional accessories

Order no.	Description	
2.142.0100	Custom Q3X thermal printer	
	Compact printer with USB interface for	
	<ul style="list-style-type: none">• 900 Touch Control• 915 KF Ti-Touch• 916 Ti-Touch• 917 Coulometer• 877 / 848 Titrino plus• 865 / 876 Dosimat plus• 91X Meter (6.2151.140)• Eco Dosimat / Titrator	
	Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.	
6.2148.050	USB/RS-232 Converter for 900 Touch Control	
	Extends the 900 Touch Control by an RS-232 interface for the connection of e.g. balance.	
6.2151.030	Cable USB A - USB B / 0.6 m	
	Connecting cable USB A - USB B	
6.2723.300	Protective cover for 900 Touch Control	
	Protective cover for 900 Touch Control	