

Accessories for 2.921.0110






Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.921.0110

Qty.	Order no.	Description
1 pcs	1.921.0110	NIRS XDS RapidContent Module
		<p>Rapid, nondestructive analyses of solid and liquid substances</p> <p>The NIRS XDS RapidContent Analyzer enables rapid, nondestructive analyses of solid and liquid formulations and substances. The sample is centered with the aid of an iris, thus enhancing the reproducibility of the analysis. A gold diffuse reflector enables the measurement of liquid samples. The Solids Module of the NIRS XDS RapidContent Analyzer extends the area of application to include practically every form of solid samples, whether it be fine powders, coarse granulates, pellets or flakes.</p>
1 pcs	6.7425.000	NIRS XDS iris
		<p>Iris for focusing mini sample vessels, reflection standard, various vials or a transflection vessel in the NIRS XDS RCA Analyzer, the NIRS XDS MasterLab Analyzer and the NIRS XDS MultiVial Analyzer</p>



Optional accessories

Order no.	Description	
6.7402.110	Einwegvials 15 mm Reflektion	
	123 verschliessbare Einwegvials aus Glas mit einem Durchmesser von 15 mm für Analysen von Feststoffen in Reflektion. Geeignet für folgende Analyser:	
	<ul style="list-style-type: none"> • NIRS DS2500 Analyzer • NIRS XDS RapidContent (Solid) Analyzer • NIRS XDS MultiVial Analyzer • NIRS XDS Masterlab Analyzer 	
6.7402.120	Einwegvials 19 mm Reflektion	
	225 verschliessbare Einwegvials aus Glas mit einem Durchmesser von 19 mm für Analysen von Feststoffen in Reflektion. Geeignet für folgende Analyser:	
	<ul style="list-style-type: none"> • NIRS DS2500 Analyzer • NIRS XDS RapidContent (Solid) Analyzer • NIRS XDS MultiVial Analyzer • NIRS XDS Masterlab Analyzer 	
6.7402.130	Disposable vials 22 mm Reflection	
	172 lockable disposable glass vials with a diameter of 22 mm for analyses of solids in reflection. Suitable for the following analyzers:	
	<ul style="list-style-type: none"> • NIRS DS2500 Analyzer • NIRS XDS RapidContent (Solid) Analyzer • NIRS XDS MultiVial Analyzer • NIRS XDS MasterLab Analyzer 	
6.7402.140	Disposable vials 28 mm Reflection	
	216 lockable disposable glass vials with a diameter of 28 mm for analyses of solids in reflection. Suitable for the following analyzers:	
	<ul style="list-style-type: none"> • NIRS DS2500 Analyzer • NIRS XDS RapidContent (Solid) Analyzer • NIRS XDS MultiVial Analyzer • NIRS XDS MasterLab Analyzer 	
6.7402.200	Stopper for 6.7402.110	
	123 stoppers for 15 mm Reflection disposable vials.	

6.7402.210 **Stopper for 6.7402.120**

225 stoppers for 19 mm Reflection disposable vials.



6.7402.220 **Stopper for 6.7402.130**

172 stoppers for 22 mm Reflection disposable vials.



6.7402.230 **Stopper for 6.7402.140**

216 stoppers for 28 mm Reflection disposable vials.

