

## Accessories for 6.01032.430

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.



### Scope of delivery 6.01032.430

Qty.	Order no.	Description
2 pcs	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optional accessories

Order no.	Description
6.01032.500	<b>Metrosep A Supp 17 Guard/4.0</b> <p data-bbox="129 472 1501 562">The Metrosep A Supp 17 Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. Thanks to the "On Column Guard System", the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column. No tools are required.</p> 
6.01032.510	<b>Metrosep A Supp 17 S-Guard/4.0</b> <p data-bbox="129 763 1501 853">The Metrosep A Supp 17 S-Guard/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is connected to the separation column simply using capillary connections. No tools are required.</p> 
6.01032.530	<b>Metrosep A Supp 17 S-Guard - 50/4.0</b> <p data-bbox="129 1055 1501 1144">The Metrosep A Supp 17 S-Guard 50/4.0 reliably protects the Metrosep A Supp 17 analytical separation columns against contamination. The guard column is simply connected to the separation column using capillary connections. No tools are required.</p> 