





Accessories for 6.03001.220

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.03001.220

Qty.	Order no.	Description	
1 pcs	6.01503.220	OMNIS 20 mL cylinder unit without accessories	
		Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Without dosing tubing and antidiffusion tip.	
1 pcs	6.1446.030	Ball stopper / B-14/(15)	
		For the movable holding of buret tips in titration vessels	
2 pcs	6.1446.040	Threaded stopper / M6	
			
1 pcs	6.1543.200	Titration tip / M6 thread	
		Titration tip with antidiffusion valve.	

2 pcs **6.1805.100** **FEP tubing / M6 / 40 cm**
With light and kink protection



1 pcs **6.2008.030** **Buret tip storage vessel for Dosing Unit**
Together with 6.2008.050 buret tip holder provides a fastening for the buret tips on 6.3032.XXX Dosing units.



1 pcs **6.2709.070** **Guiding sleeve**
For fixation of buret and dosing tips in titration heads and electrode holders



1 pcs **6.2739.000** **Wrench**
For tightening connectors



1 pcs **6.2803.010** **Grease (2 g)**
Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Optional accessories

Order no.	Description
6.05700.140	Consumable Kit for OMNIS cylinder unit

FEP M6 (40 cm) tubing connections and antidiffusion tip for OMNIS 6.03001.xxx cylinder units.

