




Accessories for 6.0331.010BR

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.0729.100	Ag/AgCl reference electrode (length 12.5 cm)	
	Silver / silver chloride reference electrode shielded with titanium shield with double junction system, installation length 10 cm. This reference electrode is well suited for: <ul style="list-style-type: none">• non-aqueous samples• differential potentiometry The SGJ 14/15 enables easy assembly and the flexible ground-joint diaphragm, which is insensitive to contamination, can be easily replaced. The reference and bridge electrolyte can be freely selected according to the application and are easy to exchange. This sensor is delivered without electrolyte filling.	
6.1241.030	Electrode shaft	
	Electrode shaft for separate metal pin electrodes	
6.2106.020	Strand / 1 m / 2 x B	
	Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2310.010	Electrolyte KNO₃-1mol/L 250 mL	
	Electrolyte solution KNO ₃ 1 M (reference electrolyte for combined silver electrodes)	