






Accessories for 6.05609.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05609.000

Qty.	Order no.	Description	
1 pcs	6.01405.010	Titration vessel lid for volumetric Karl Fischer titration Titration vessel lid for volumetric Karl Fischer titration with an OMNIS titrator.	
1 pcs	6.01406.010	Adsorber tube with lid Adsorber tube with lid for OMNIS Karl Fischer titrator, volumetric	
1 pcs	6.01406.220	20 – 90 mL titration vessel 20 – 90 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.	
1 pcs	6.01406.250	50 – 150 mL titration vessel 50 – 150 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.	
1 pcs	6.01408.010	Set of seals for Karl Fischer titration cell Complete set of seals, suitable for 6.05609.000 and 6.05609.010.	

1 pcs 6.01543.120 Antidiffusion tip

Titration tip with antidiffusion valve, suitable for OMNIS Karl Fischer titrator.



1 pcs 6.02104.040 Electrode cable plug-in head G (pol.) / plug P, 0.55 m

Cable for connecting an analog polarizable electrode (plug-in head G) to an analog measuring module from OMNIS.



3 pcs 6.02709.010 Stopper

Plastic stoppers, suitable for 2.1001.4220 OMNIS Karl Fischer Titrator.



1 pcs 6.02709.020 Septum stopper

Plastic septum stopper, suitable for 6014050X0 titration vessel lids.



2 pcs 6.02709.030 Screw nipple

Plastic screw nipple, suitable for 6014050X0 titration vessel lids.



1 pcs 6.02801.000 Cotton for adsorber tube

Wad of cotton (4 pieces), suitable for OMNIS 601406010 adsorber tube.



1 pcs **6.0338.100** **Double Pt-wire electrode for volumetry**

Indicator electrode used for volumetric Karl Fischer titration.



2 pcs **6.1448.010** **Septum 12 mm, 5 pieces**

Set of 5 items.



2 pcs **6.1903.020** **Stirring bar / 16 mm**

Stirring Bar with magnetic core, PTFE coated, length 16 mm



2 pcs **6.1903.030** **Stirring bar / 25 mm**

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs **6.2412.000** **Weighing boats**

Weighing boat made of glass for the 6.1414.030 titration vessel lid or the 6.1455.31X, 6.1464.32X, 6.1465.320 titration vessel. For KF water determinations. Including protective tube.



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



2 pcs 6.2816.020 Syringe 10 mL with Luer connection

For various applications in IC and VA



2 pcs 6.2816.030 Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.

