


Accessories for 6.05700.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.05700.010

Qty.	Order no.	Description
2 pcs	6.01801.020	Pharm-A-line pump tubing to the peristaltic pump module
		The pump tubing comes pre-equipped with adapters that permit direct installation in the peristaltic pump of the pump module. The Pharm-A-line pump tubing is used for the transport of both aqueous and nonaqueous solutions and must be operated in connection with the white press clamp in the pump module.
		
2 pcs	6.02703.000	Press clamp, white, for "Peristaltic" pump module
		The press clamp is used to fix and secure the pump tubing in the "Peristaltic" pump module.
		