



Accessories for 6.0750.100

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.0750.100

Qty.	Order no.	Description	
1 pcs	6.1236.050	Sleeve with SGJ 14/12 mm Sleeve with SGJ 14/12 mm, polyethylene.	
1 pcs	6.2008.040	Storage vessel Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units.	

Optional accessories

Order no.	Description	
6.2103.140	Adapter black 2 mm plug / B socket 4 mm For connecting plug B (4 mm) to 2 mm socket.	
6.2106.020	Strand / 1 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2106.060	Strand / 2 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	

6.2310.010 Electrolyte KNO₃-1mol/L 250 mL

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)



6.2320.000 TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), c(TEABr) = 0.4 mol/L

