



## Accessories for 6.0915.100

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.0915.100

| Qty.  | Order no.  | Description  |  |
|-------|------------|--|--|
| 1 pcs | 6.1236.050 | <b>Sleeve with SGJ 14/12 mm</b><br>Sleeve with SGJ 14/12 mm, polyethylene.   |   |
| 1 pcs | 6.2008.040 | <b>Storage vessel</b><br><br>Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units. |  |

## Optional accessories

| Order no. | Description |
|-----------|-------------|
|-----------|-------------|

**6.2301.060 KCl ion standard 250 mL**

KCl ion standard ( $c(\text{ion}) = 0.1000 \pm 0.0005 \text{ mol/L}$ ), can also be used as standard solution for conductometry  $12.87 \text{ mS/cm}$  ( $25^\circ\text{C}$ ).



**6.2324.010 Conductivity standard 100  $\mu\text{S/cm}$  250 mL**

Conductivity standard for calibration of conductivity measuring cells with cell constant =  $0.1/\text{cm}$ .



**6.2324.110 Conductivity standard 100  $\mu\text{S/cm}$ , 5 x 30 mL**

Conductivity standard for calibration of conductivity measuring cells with cell constant =  $0.1/\text{cm}$ .



**6.2763.000 Flow-through cell for the 6.0915.100/6.0920.100**

Flow-through cell for the 6.0915.100/6.0920.100 5-ring conductivity measuring cell

