

## Accessories for 6.1010.410

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description
-----------	-------------

<b>6.1010.450</b>	<b>Metrosep C 3 Guard/4.0</b>
-------------------	-------------------------------

The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be prolonged considerably. The Metrosep C 3 Guard/4.0 also functions according to the "On Column Guard System" and is screwed directly onto the separation column with nearly no dead volume.



---

<b>6.1010.460</b>	<b>Metrosep C 3 S-Guard/4.0</b>
-------------------	---------------------------------

The Metrosep C 3 S-Guard/4.0 contains the C-3 column material and is used to protect Metrosep C-3 cation columns. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be prolonged considerably.

