

Accessories for 6.1020.220

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.




Scope of delivery 6.1020.220

Qty.	Order no.	Description
2 pcs	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.1020.600	Metrosep A Supp 10 Guard/2.0	
	<p>The Metrosep A Supp 10 Guard/2.0 protects the analytical A Supp 10 microbore columns effectively from contamination. Thanks to its «On Column Guard System» the Guard Column is characterized by its easy handling. It is simply screwed onto the analytical column; no tools are required.</p>	
6.1020.610	Metrosep A Supp 10 S-Guard/2.0	
	<p>The Metrosep A Supp 10 S-Guard/2.0 protects the analytical A Supp 10 microbore columns effectively from contamination.</p>	