







Accessories for 6.1459.310

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.01403.060	Titration head 3xNS14 / 4x6.4 mm for 200 mL disposable beaker (P&P)	
	Titration head with 3 NS14 and 4x 6.4 mm openings for using 200 mL PP disposable beakers (6.1459.310) or 150 mL beaker glasses (high shape, without spout) in Pick&Place modules of the OMNIS Sample Robot Pick&Place with two retracted dosing tips. Prepared for being equipped with rinsing and aspiration tubing.	
6.01404.020	Beaker adapter for 200 mL sample beaker	
	Beaker adapter for centering the 200 mL sample beaker in the slide of the Pick&Place modules.	
6.02041.020	OMNIS sample rack 9 x 200 mL	
	OMNIS 9 x 200 mL sample rack for use with disposable PP sample beaker.	
6.02601.020	Gripper fingers 50 - 72 mm	
	Gripper fingers for OMNIS Sample Robot Pick&Place for gripping sample beakers with an outer diameter of 50 to 72 mm.	
6.1458.010	Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings	
	With 6xSGJ 14 and 3x SGJ 9 apertures.	
6.1909.010	Stirring propeller / 96 mm	
	Propeller stirrer, length from lower edge of SGJ: 96 mm. For use with beakers with 722, 802 Rod Stirrer	

6.1909.020 Stirring propeller / 104 mm

Propeller stirrer, length from lower edge of SGJ: 104 mm. For use with beakers with 722, 802 Rod Stirrer



6.1909.040 Stirring propeller / 102 mm

Propeller stirrer, length from lower edge of SGJ: 102 mm. For titration vessels with SGJ 14/15 and 722, 802 Rod Stirrer



6.2041.370 Sample rack 14 x 200 mL

Sample rack for 14 x 200 mL samples, one-row



6.2041.830 Sample rack 28 x 200 mL

Sample rack for 28 x 200 mL samples (6.1459.310), multi-row

