

Accessories for 6.1826.060

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description	
6.2744.030	Coupling olive/UNF 10/32, 4x	
	Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump	
6.2744.110	Pump tubing - inline filter PEEK adapter	
	Connector between pump tubing and inline filter	
6.2744.160	Nozzle/UNF 10/32 coupling with security device	
	Security device for connecting pump tubing to a nozzle.	
6.2744.180	Pump tubing connector with security lock and filter	
	For the connection of a pump tubing and a capillary with built-in filter	