


Accessories for 6.2125.090

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

| Order no. | Description | |
|--|---------------------------------------|---|
| 6.2125.120 | Adapter sub D 25-pin/5 sockets |  |
| Distributes the «Remote lines to 5 x 25-pin sockets. Used with pH Meters, Titrinos, KF Coulometer, Sample Changers, IC systems, etc. | | |
