

## Accessories for 6.2729.110

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.


### Scope of delivery 6.2729.110

Qty.	Order no.	Description
3 pcs	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



## Optional accessories

Order no.	Description	
<b>6.2744.000</b>	<b>Pressure screw PVDF</b>	
With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)		

---