





Accessories for 6.2824.040

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2824.040

Qty.	Order no.	Description	
2 pcs	6.2741.020	Piston seal (orange) For all standard pump heads	
1 pcs	6.2824.010	Outlet valve (steel) Spare part for IC Pumps (steel pump head)	
1 pcs	6.2824.020	Inlet valve (steel) Spare part for IC Pumps (steel pump head)	

Optional accessories

Order no.	Description	
6.2824.050	Spring for main piston For pump heads 6.2824.100	
6.2824.060	Spring for all pump pistons For all pump heads	
6.2824.070	Zirconium oxide piston Pump piston for all standard pump heads	