





Accessories for 6.2841.130

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.2841.130




Qty.	Order no.	Description	
1 pcs	6.1801.120	PVC tubing / 4 mm / 6 mm / 2 m	
1 pcs	6.1808.050	Tubing adapter olive / M8 outer 1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.	
1 pcs	6.2621.030	Hexagon key 4 mm	
1 pcs	6.2621.140	Hexagon key 2.5 mm	

2 pcs **6.2744.060** **Threaded stopper**

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Optional accessories

Order no.	Description	
6.2762.100	Mixing vessel for Liquid Handling Station	
	Replacement mixing vessel for the combined rinsing and dilution stations.	
6.2762.110	Cover for the 6.2762.100 mixing vessel	
	Replacement cover for the mixing vessel.	
6.2762.200	Transparent cover for the 6.2762.100 mixing vessel, PMMA	
	Transparent cover for the 6.2762.100 mixing vessel, for visual inspection in the mixing vessel.	