





## Accessories for 6.5000.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5000.000

Qty.	Order no.	Description	
1 pcs	6.1602.210	<b>Attachment 5 x UNF 10/32 for «Waste Collector»</b> Attachment for holding up to 5 waste capillaries.	
1 pcs	6.1807.010	<b>Y-connector for tubing 6-9 mm i.d.</b> Connector for waste tubing	
1 pcs	6.1815.010	<b>Spiral band 0.5 m</b> For holding together different cables or tubing.	
2 pcs	6.1816.020	<b>Silicone tubing 6 mm i.d., 1 m</b> For waste tubing.	

---

**1 pcs**    **6.1816.080**    **Silicone tubing 6 mm i.d. / 2 m**

Drainage tubing to the "Waste Collector".



---

**1 pcs**    **6.2057.220**    **Holder to Waste Collector**

Holder for fastening the Waste Collector to the IC instrument.



---

**1 pcs**    **6.2151.020**    **Cable USB A - USB B 1.8 m**  
USB connecting cable



---

**1 pcs**    **6.2251.000**    **Colored sleeves for capillaries**  
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



---

**1 pcs**    **6.2617.010**    **Tool for piston seal**  
For removing and assembling the piston seal for all standard pump heads



---

**2 pcs**     **6.2621.000**     **Adjustable wrench**  
Max. opening: 20 mm. For IC instruments



---

**1 pcs**     **6.2621.030**     **Hexagon key 4 mm**



---

**1 pcs**     **6.2621.050**     **1/4 in. wrench**  
For 1/4 in. screws. For IC instruments



---

**1 pcs**     **6.2621.080**     **Capillary cutter**  
For plastic capillaries. Used with IC instruments



---

**1 pcs**      **6.2621.100**      **3 mm hex key for IC Sample Processors**



---

**2 pcs**      **6.2739.000**      **Wrench**  
For tightening connectors



---

**1 pcs**      **6.2743.080**      **Stopper for overflow (ProfIC), 5 pieces**  
For Professional IC instruments



---

**1 pcs**      **6.2765.000**      **Container with GL 45 drainage to the "Waste Collector"**

Container with drainage to the "Waste Collector".

