

Accessories for 6.5327.020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5327.020

Qty.	Order no.	Description	
1 pcs	6.0343.000	Separate Pt Rod Electrode Platinum auxiliary electrode for voltammetry and CVS.	
1 pcs	6.0728.030	LL-Ag/AgCl reference electrode for CVS Inner reference system with electrolyte c(KCl) = 3 mol/L. For use with electrolyte vessel made of plastic 6.1245.010.	
1 pcs	6.1204.210	Driving axle for rotating disk electrode (RDE) Driving axle for rotating disk electrode (RDE) with titanium axle for VA measuring stands, M3 thread.	
1 pcs	6.1204.610	2 mm platinum electrode tip for CVS Metal electrode tip made of platinum (Pt), diameter of electrode disk 2 mm ±0.05 mm, polished surface, shaft made of glass, M3 thread connection. Concentricity error: < 0.3 mm. Temperature range: 0 - 50 °C.	
1 pcs	6.1245.010	Electrolyte vessel made of plastic Electrolyte vessel made of plastic with ceramic diaphragm for double junction reference electrode.	

1 pcs **6.1415.210** **Titration vessel / 10-90 mL**



1 pcs **6.1415.250** **Titration vessel / 50-150 mL**



1 pcs **6.2310.010** **Electrolyte KNO₃-1mol/L 250 mL**
Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)



1 pcs **6.2709.040** **Stopper for MME opening on VA Stand**
For closing gas tubing of the Multi-Mode Electrode in VA Stands if not in use. E.301.0004 O-rings included.

