





Accessories for 6.5335.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5335.000

Qty.	Order no.	Description	
2 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32	
		For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1608.030	Round glass bottle / 1000 mL / GL 45	
1 pcs	6.1618.050	Thread adapter 40 mm to GL 45 Adapter for bottles with 40 mm thread to GL45.	
1 pcs	6.1621.000	Container 10 L As rinsing or waste container in automated systems.	

1 pcs 6.1826.360 Pump tubing LFL (white/white), 3 stoppers
For Sample Processors



2 pcs 6.1826.390 Pump tubing LFL (yellow/yellow), 3 stoppers
For the sample solution in Inline Ultrafiltration.



3 pcs 6.1826.420 PharMed® pump tubing (orange/yellow), 3-stopper

For MSM regeneration with peristaltic pump.



4 pcs 6.2744.010 Pressure screw 5x
With UNF 10/32 connection. For the connection of PEEK capillaries



2 pcs 6.2744.030 Coupling olive/UNF 10/32, 4x
Connection between pressure screw and pump tubing. 4 pieces. For IC instruments with peristaltic pump



6 pcs 6.2744.160 Nozzle/UNF 10/32 coupling with security device
Security device for connecting pump tubing to a nozzle.

