





## Accessories for 6.5610.010

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 6.5610.010

Qty.	Order no.	Description	
1 pcs	6.0340.000	<b>Double Pt-wire electrode for sample changer (fixed cable)</b>  Indicator electrode (diameter = 5.3 mm, fixed cable 2 m) used for automated Karl Fischer titration.	
24 pcs	6.1432.210	<b>Sample beaker, clear glass, 75 mL, 1 piece</b>  Sample beaker, clear glass, 75 mL, suitable for sample racks 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece	
1 pcs	6.1458.030	<b>KF titration head / 4 x M10 connectors</b>  With 4 M10 screw nipples, including stoppers.	
1 pcs	6.1805.060	<b>FEP tubing / M6 / 60 cm</b>  With light and kink protection	

---

**1 pcs**    **6.1808.000**    **Connection bushing M6**



---

**1 pcs**    **6.1820.020**    **M10/M6 adapter**

Connector with M10 thread for 6.1812.XXX tubing and M6 thread for 6.1805.XXX tubing.



---

**1 pcs**    **6.1821.000**    **Aspiration tube**

For KF Sampler



---

**24 pcs**    **6.1903.010**    **Stirring bar / 12 mm**

Stirring Bar with magnetic core, PTFE coated, length 12 mm



---

**24 pcs**    **6.2037.040**    **Foil holder**

Foil holder for 75 mL sample beaker for Karl Fischer for KF Sample Changer



---

**1 pcs**     **6.2041.340**     **Sample rack 24 x 75 mL**  
Sample rack for 24 x 75 mL samples, one-row.



---

**1 pcs**     **6.2820.000**     **Al foil**  
1000 circular sheets of 80 mm diameter. Thickness: 0.010 mm. Used together with 6.2037.030 Foil holder for 6.1432.210 Sample beaker. For 674 KF Sample Changer

