

Accessories for 6.5613.000

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5613.000

Qty.	Order no.	Description	
1 pcs	6.1403.040	KF Adsorber tube Used with KF equipment. Including cover and O-ring.	
1 pcs	6.1414.060	Micro titration vessel lid ECO	
1 pcs	6.1415.150	Titration vessel / 5-70 mL	
2 pcs	6.1446.040	Threaded stopper / M6	

1 pcs **6.1448.010** **Septum 12 mm, 5 pieces**
Set of 5 items.



1 pcs **6.1543.120** **Aspiration tip / M8 thread**
Aspiration tip if 6.1805.200 siphoning tubing is used



2 pcs **6.1805.010** **FEP tubing / M6 / 13 cm**
With light and kink protection.



1 pcs **6.1805.200** **PTFE tubing / M8 / 0.5 m**



1 pcs **6.1808.100** **Distributor 3 x M6 / 1 x M8**
For 6.5613.000 ECO titration equipment



1 pcs **6.1903.020** **Stirring bar / 16 mm**
Stirring Bar with magnetic core, PTFE coated, length 16 mm



2 pcs **6.2709.090** **Stopper**
PTFE. Used with 6.5613.000 ECO titration equipment



1 pcs **6.2730.020** **Septum stopper**
With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.



5 pcs **6.2730.030** **Stopper**
Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



2 pcs **6.2730.060** **Screw nipple**
PCTFE. With M16 thread. Used with 6.5613.000 ECO titration equipment



2 pcs **6.2730.070** **Screw nipple**
PCTFE. With M16 thread. Used with 6.5613.000 ECO titration equipment



1 pcs **6.2739.000** **Wrench**
For tightening connectors



2 pcs **6.2740.030** **Rinsing nozzle**
The rinsing solution is applied as a fine jet for the optimal removal of layers on the electrodes and titration accessories. For 6.5613.000 ECO titration equipment.



Optional accessories

Order no.	Description
6.2730.080 PCTFE. Used with 6.5613.000 ECO titration equipment with 722 Stirrer	Screw nipple

