

Accessories for 6.7401.320

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Optional accessories

Order no.	Description
-----------	-------------

6.7403.110 **NIRS XDS spacer for 0.5 mm cuvettes (6.7401.110)**

Spacer for fixing open quartz cuvettes with 0.5 mm total pathlength in an NIRS XDS RLA Analyzer.



6.7403.180 **NIRS Spacer for 12.5 mm cuvettes**

The NIRS Spacer for 12.5 mm cuvette set ensures spacers for an optimum position of the cuvette in the sample holder.

