

Accessories for 6.7402.030

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.7402.030

Qty.	Order no.	Description
1 pcs	6.7402.020	NIRS disposable backs for 6.7402.030, 100 pcs.

The disposable backs are used for the secure sealing of the mini sample cups. The sample cups are used in the following:

- NIRS DS2500 Analyzer (order number: 2.922.0010)
- NIRS XDS MasterLab Analyzer (order number: 2.921.1310)
- NIRS XDS RCA Analyzer (order number: 2.921.1110)
- NIRS XDS MultiVial Analyzer (order number: 2.921.1210)

The package is comprised of 100 disposable backs.

