





## Accessories for FLWCL-SC

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description	
<b>110</b>	<b>Screen-Printed Carbon Electrode</b>	
	Screen-Printed Carbon Electrode (Aux.:C; Ref.:Ag). Suitable for working with microvolumes, for decentralized assays or to develop specific sensors	
<b>150</b>	<b>Screen-Printed Carbon Electrode (Aux.:Pt; Ref.:Ag)</b>	
	Screen-Printed Carbon Electrode (Aux.:Pt; Ref.:Ag). Suitable for working with microvolumes, for decentralized assays or to develop specific sensors.	
<b>CAC</b>	<b>Cable connector for Screen-Printed Electrodes</b>	
	Connects DropSens Screen-Printed Electrodes to any kind of potentiostat. Recommended when working in batch or flow analysis.	
<b>CAST</b>	<b>µStat Cable connector for Screen-Printed Electrodes</b>	
	Connects any DropSens equipment with DropSens Screen-Printed Electrodes	