






Accessories for 2.945.8020

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.945.8020

Qty.	Order no.	Description	
1 pcs	1.891.0010	891 Professional Analog Out	
		The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.	
1 pcs	1.942.0300	942 Extension Module Vario ChS/PP	
		The 942 Extension Module Vario ChS/PP is an extension module for the Professional IC Vario line of instruments. It permits the integration of a complete chemical suppression with a peristaltic pump for suppressor regeneration or Inline Neutralization or Inline Cation Removal , respectively, in a Professional IC Vario system.	
		Typical applications:	
		<ul style="list-style-type: none"> • Chemical suppression for 2nd analysis channel • Inline Neutralization/Cation Removal 	
1 pcs	1.945.0010	945 Professional Detector Vario – Conductivity	
		Intelligent stand-alone detector equipped with a high-performance IC Conductivity Detector. For use with intelligent IC instruments or as independent conductivity detector.	
1 pcs	6.1602.150	Bottle attachment / GL 45 - 3 x UNF 10/32	
		For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1608.020	Glass bottle / 1000 mL / GL 45	
		Bottle for auxiliary solutions.	

1 pcs **6.1803.030** **PTFE capillary 0.5 mm i.d. / 3 m**
Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



2 pcs **6.1826.420** **PharMed® pump tubing (orange/yellow), 3-stopper**

For MSM regeneration with peristaltic pump.



1 pcs **6.1831.030** **PEEK capillary 0.75 mm i.d. / 3 m**

For all IC components.



1 pcs **6.2061.100** **Bottle holder for Professional IC instruments**

Instrument cover and bottle holder for Professional IC instruments.



1 pcs **6.2061.110** **Tray with sensor for Professional IC Instruments**



1 pcs **6.2128.200** **Cable 891-Analog**
Analog cable to connect to external AD converters.



1 pcs **6.2141.380** **Remote cable Professional IC - MS-Detector (Agilent)**
Cable for synchronization of MagIC Net™ and Chemstation (Agilent).



1 pcs **6.2148.010** **Remote Box MSB**
Additional remote interface for the connection of devices that can be controlled via remote lines.
With permanently attached cable.



1 pcs **6.2151.020** **Cable USB A - USB B 1.8 m**
USB connecting cable



2 pcs **6.2156.070** **Cable Extension Module - Professional IC, 1 m**
Cable connecting the extension module with the Professional IC instrument



1 pcs **6.2251.000** **Colored sleeves for capillaries**
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 pcs **6.2744.014** **Pressure screw 2x**
 With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**
 Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



2 pcs **6.2744.180** **Pump tubing connector with security lock and filter**
 For the connection of a pump tubing and a capillary with built-in filter



1 pcs **6.2842.000** **MSM-HC Rotor A**
 Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)



1 pcs **6.5000.030** **Accessory kit: Vario/Flex ChS**

Set of accessories for chemical suppression (ChS) for Professional IC Vario and Compact IC Flex instruments. Without Suppressor Rotor.



1 pcs 6.6059.322 MagIC Net 3.2 Professional CD: 1 license



Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.


Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.

MagIC Net Professional complies fully with FDA Regulation 21 CFR Part 11 as well as GLP.

MagIC Net is available in 16 dialog languages.

1 license.

Optional accessories

Order no.	Description
6.6059.323	MagIC Net 3.2 Multi CD: 3 licenses
<p data-bbox="1268 414 1516 504"></p> <p data-bbox="130 481 1252 582">Professional PC program for controlling all of the intelligent Professional IC systems, Compact IC instruments and their peripherals, such as various Autosamplers, 800 Dosinos, 771 Compact Interface, etc. The software permits the inspection, data acquisition, data evaluation, data monitoring and report generation of ion chromatographic analyses.</p> <p data-bbox="130 604 1252 705">Graphical user interface for routine operations, extensive database programs, method development, configuration and manual system control; very flexible user administration, efficient database operations, extensive data export functions, individually configurable report generator, control and monitoring of all system components and the chromatography results.</p> <p data-bbox="130 728 997 761">MagIC Net Multi complies fully with FDA Regulation 21 CFR Part 11 as well as with GLP.</p> <p data-bbox="130 772 598 806">MagIC Net is available in 16 dialog languages.</p> <p data-bbox="130 817 502 851">Client/server version with 3 licenses.</p>	
