


Accessories for 2.136.0800

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.136.0800

Qty.	Order no.	Description
1 pcs	1.136.0800	MMS 5000 Autosampler
		MMS 5000 Autosampler from Analytik Jena for use with the Metrohm Combustion IC for fully automatic analysis of liquid and solid samples. In order to match the modular Multi-Matrix sampler with the correct sample type, either the liquid kit (6.7303.000) or the solid kit (6.7302.000) must be used.



Optional accessories

Order no.	Description	
6.7301.150	Substrate made of quartz for liquid samples	
	Substrate for liquid sample injections, 3 pcs.	
6.7302.000	Kit for solid samples for MMS 5000, 35 positions	
	Kit for solid samples for MMS 5000 with rack for sample boats, grippers and 10 quartz glass boats.	
6.7302.100	Quartz boat 40 x 9 mm	
	Quartz glass boat for sample injection in Combustion Module, 10 pcs.	
6.7303.000	Kit for liquid samples for MMS 5000, 112 positions	
	Kit for liquid samples for MMS 5000 with rack for sample vessels, 100 µL syringe with replaceable needle and 100 snap-ring bottles.	
6.7303.100	Snap-ring bottle ND11, 1.5 mL	
	For liquid holder with MMS 5000, 100 pcs.	
6.7303.110	Septum cap ND11	
	For liquid holder with MMS 5000, 100 pcs.	
6.7303.200	50 µL syringe with replaceable needle	
	Replacement syringe for MMS 5000 liquid setup	

6.7303.210 100 μ L syringe with replaceable needle

Replacement syringe for MMS 5000 liquid setup



6.7303.220 Needles for 50 μ L/100 μ L syringes, 3 pcs.

Replacement needle for syringe for MMS 5000 liquid setup

