






Accessories for 6.5330.170

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 6.5330.170

Qty.	Order no.	Description	
1 pcs	6.1014.200	Metrosep I Trap 1 - 100/4.0	
		The Metrosep I Trap 1 - 100/4.0 column is used to eliminate ionic, i.e. cationic and anionic, contaminants from aqueous solutions. Its use is especially recommended for the purification of the transfer water in combination with "MISP" or "Metrohm Inline Sample Preparation." Using the Metrosep I Trap 1 - 100/4.0 column can significantly reduce the influence of the transfer water on the system blank.	
1 pcs	6.1602.160	Bottle attachment GL 45 for eluents and reagents	
		For eluent, reagent, and rinsing bottles, with connectors for the adsorber tube and the aspiration tubing.	
1 pcs	6.1608.070	Eluent bottle / 2 L / GL 45	
1 pcs	6.1618.020	Thread adapter S 40 to GL 45	
		Adapter for bottles with S40 thread (Merck...) to GL45.	
1 pcs	6.1624.000	Adapter SGJ 14 for 6.1619.XXX Adsorber tube	
		Adapter SGJ 14. Part of 6.1619.010 for Dilution Sample Processor.	

1 pcs **6.1805.120** **FEP tubing / M6 / 100 cm**
With light and kink protection



1 pcs **6.1808.090** **Thread adapter M8 outer / M6 inner**
Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



1 pcs **6.1819.110** **FEP aspiration tubing to canister**
Aspiration tubing to canister attachment with M6 connector. For the production of eluents.



1 pcs **6.1825.290** **PEEK sample loop 250 µL**
Sample loop for «Metrohm intelligent Partial loop Technique» (MiPT) and for PILS (particle into liquid sampler). With 2 Peek pressure screws.



1 pcs **6.1831.100** **PEEK capillary 0.25 mm i.d., 1 m**



1 pcs **6.1841.000** **PEEK transfer capillary 2 mL, 5 m**



1 pcs **6.2057.210** **Holder for Dosino to IC instruments**
Mounting part for Dosinos at the IC instrument housing.



1 pcs **6.2624.100** **Pick-up needle MF**
Coated steel needle for metalfree pick-up sampling as well as for Eppendorf sample vials.



1 pcs **6.2744.010** **Pressure screw 5x**
 With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.080** **M6 thread / UNF 10/32 adapter**
 With M6 outer thread. For IC instruments.



1 pcs **6.2833.010** **Needle/tube holder 1/16 in.**
 For IC Sample Processors. For 1/16 in. tubes



1 pcs **6.2841.120** **Liquid Handling Station, left**

Combined rinsing and dilution station for sample changers with Swing Head. For assembly to the left of the tower.



1 pcs **6.3032.120** **Dosing Unit 2 mL**

Dosing unit with integrated data chip with 2 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 glass thread. FEP tubing connection, anti-diffusion buret tip.



1 pcs **6.9920.231** **Retaining plate for sample rack 159 x 2 mL + 3 x 300 mL**

Plate for deflecting the vessel cover on the Professional IC Sample Processor. For 2 mL sample vessels.

