





Accessories for 2.942.0020

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.942.0020

Qty.	Order no.	Description	
1 pcs	1.942.0020	942 Extension Module Vario Prep 2	
		<p>The 942 Extension Module Vario Prep 2 is an extension module for the Professional IC Vario line of instruments. It multiplies the possibilities of Inline Sample Preparation and Liquid Handling in Professional IC Vario systems.</p> <p>Typical applications</p> <ul style="list-style-type: none"> • Inline Dialysis in two-channel Vario systems • Inline Dialysis in systems without sample changers • Column switching and cutting techniques 	
1 pcs	6.1803.020	PTFE capillary 0.97 mm i.d. / 5 m	
		For all IC instruments	
1 pcs	6.1825.230	PEEK sample loop 10 µL	
		For injection valve, with 2 PEEK pressure screws	
2 pcs	6.1826.340	Pump tubing LFL (black/black), 3 stoppers	
		For the sample solution in inline dialysis.	

1 pcs **6.1831.060** **PEEK capillary 0.5 mm i.d., 1 m**
For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit



1 pcs **6.1831.100** **PEEK capillary 0.25 mm i.d., 1 m**



1 pcs **6.2156.060** **Cable Extension Module - Professional IC, 40 cm**
Cable connecting the extension module with a Professional IC instrument



1 pcs **6.2251.000** **Colored sleeves for capillaries**
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 pcs **6.2621.100** **3 mm hex key for IC Sample Processors**



1 pcs **6.2744.010** **Pressure screw 5x**
With UNF 10/32 connection. For the connection of PEEK capillaries



1 pcs **6.2744.034** **Coupling olive/UNF 10/32, 2x**
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.



2 pcs **6.2744.070** **Pressure screw short**
Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries



2 pcs **6.2744.160** **Nozzle/UNF 10/32 coupling with security device**
Security device for connecting pump tubing to a nozzle.



Optional accessories

Order no.	Description	
6.2057.120	Ultrafiltration cell holder	
	Holder for the ultrafiltration cell	
6.2057.130	Dialysis cell holder	
	Holder for the dialysis cell	
6.2103.170	Adapter to Leak Sensor at Professional IC	
	Adapter to connect two leak sensors to a Professional IC instrument	
6.2156.070	Cable Extension Module - Professional IC, 1 m	
	Cable connecting the extension module with the Professional IC instrument	
6.5332.000	IC tubing kit for second extension module	