








## Accessories for 2.942.0300

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

### Scope of delivery 2.942.0300

Qty.	Order no.	Description	
1 pcs	1.942.0300	<b>942 Extension Module Vario ChS/PP</b>	
		<p>The 942 Extension Module Vario ChS/PP is an extension module for the Professional IC Vario line of instruments. It permits the integration of a complete <b>chemical suppression</b> with a <b>peristaltic pump</b> for suppressor regeneration or <b>Inline Neutralization</b> or <b>Inline Cation Removal</b>, respectively, in a Professional IC Vario system.</p> <p>Typical applications:</p> <ul style="list-style-type: none"> <li>• Chemical suppression for 2nd analysis channel</li> <li>• Inline Neutralization/Cation Removal</li> </ul>	
1 pcs	6.1602.150	<b>Bottle attachment / GL 45 - 3 x UNF 10/32</b>	
		For the connection of 1/16 in. capillaries. Used in conjunction with MSM auxiliary solutions and in Metrohm Inline Sample Preparation.	
1 pcs	6.1608.020	<b>Glass bottle / 1000 mL / GL 45</b>	
		Bottle for auxiliary solutions.	
1 pcs	6.1803.030	<b>PTFE capillary 0.5 mm i.d. / 3 m</b>	
		Capillary for inline dialysis, for Dialysis Unit, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit	

2 pcs	6.1826.420	<b>PharMed® pump tubing (orange/yellow), 3-stopper</b>	
For MSM regeneration with peristaltic pump.			
1 pcs	6.2156.070	<b>Cable Extension Module - Professional IC, 1 m</b>	
Cable connecting the extension module with the Professional IC instrument			
1 pcs	6.2251.000	<b>Colored sleeves for capillaries</b>	
Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.			
1 pcs	6.2744.034	<b>Coupling olive/UNF 10/32, 2x</b>	
Connection between pressure screw and pump tubing. 2 pieces. For IC instruments with peristaltic pump.			
2 pcs	6.2744.180	<b>Pump tubing connector with security lock and filter</b>	
For the connection of a pump tubing and a capillary with built-in filter			
1 pcs	6.5000.030	<b>Accessory kit: Vario/Flex ChS</b>	
Set of accessories for chemical suppression (ChS) for Professional IC Vario and Compact IC Flex instruments. Without Suppressor Rotor.			

## Optional accessories

Order no.	Description	
6.2103.170	<b>Adapter to Leak Sensor at Professional IC</b> Adapter to connect two leak sensors to a Professional IC instrument	
6.2832.000	<b>MSM Rotor A</b> Suppressor rotor for all IC instruments with MSM (Metrohm Suppressor Module)	
6.2835.000	<b>SPM Rotor A</b> Rotor for inline sample preparation with cation exchange	
6.2842.000	<b>MSM-HC Rotor A</b> Suppressor rotor for all IC instruments with MSM-HC (Metrohm Suppressor Module with high capacity)	
6.2842.020	<b>Adapter sleeve for Suppressor Vario</b> The use of the adapter sleeve makes it possible to install both the Standard MSM and the MSM-LC rotor in the MSM-HC housing.	

---

**6.2844.000**

**MSM-LC Rotor A**

Suppressor rotor for use with 2 mm columns.

