










Accessories for 2.899.0110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.899.0110

Qty.	Order no.	Description	
1 pcs	1.899.0010	899 Coulometer	
1 pcs	6.0341.100	Double Pt-wire electrode for coulometry Indicator electrode used for coulometric Karl Fischer titration.	
1 pcs	6.0345.100	Generator electrode without diaphragm Generator electrode without diaphragm for coulometric Karl Fischer titrations, SGJ 29/22. This generator electrode is well suited for most applications where the sample is highly soluble in alcohol.	
1 pcs	6.1403.030	KF adsorber tube for coulometer cell For KF coulometers.	

1 pcs	6.1437.000	SGJ stopper SGJ 14	
<hr/>			
2 pcs	6.1448.020	Septum 16 mm, 5 pieces	
	Set of 5 pcs.		
<hr/>			
1 pcs	6.1464.320	KF titration vessel / 80-250 mL / coulometric	
<hr/>			
1 pcs	6.1903.030	Stirring bar / 25 mm	
		Stirring bar with magnetic core, PTFE covering, length 25 mm.	
<hr/>			
1 pcs	6.2047.040	Titration vessel holder for 899/917	
		Titration vessel holder for coulometric cells for 899 Coulometer and 917 Ti-Touch	

1 pcs **6.2103.130** **Adapter red 2 mm plug / 4 mm socket**
For connecting plug B (4 mm) to 2 mm socket.



1 pcs **6.2103.140** **Adapter black 2 mm plug / B socket 4 mm**
For connecting plug B (4 mm) to 2 mm socket.



1 pcs **6.2104.020** **Electrode cable / 1 m / F**
For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).



1 pcs **6.2104.120** **Electrode cable / 1 m / H**
With plug H. For connecting generator electrode - KF Coulometer.



1 pcs **6.2147.000** **Numerical USB keypad**
Numerical USB keypad for Titrino plus, Dosimat plus, 862 Compact Titrosampler, 869 Compact Sample Changer, 863 Compact Autosampler and 899 Coulometer. With two integrated USB connectors.



1 pcs **6.2151.100** **Adapter USB MINI (OTG) - USB A**

For connecting USB instruments.



1 pcs **6.2164.010** **Power supply unit 100 - 240 V/24 V DC**

Power supply unit for 915 KF Ti-Touch and 916 Ti-Touch



1 pcs **6.2701.040** **GL 18 screw cap with hole**

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



3 pcs **6.2713.000** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 14.



1 pcs **6.2713.010** **Sleeve**

Sleeve for the fat-free ground-joint connection with the standard ground-joint SGJ 29.



1 pcs **6.2713.020** **Sleeve**

For the fat-free ground-joint connection with the standard ground-joint SGJ 19.



1 pcs 6.2723.310 Protective cover for Titrino plus

Protective cover for display and keypad for all instruments of the Titrino plus family.

Suitable for: 848 Titrino plus, 862 Compact Titrosampler, 863 Compact Autosampler, 865 Dosimat plus, 869 Compact Sample Changer, 870 KF Titrino plus, 876 Dosimat plus, 877 Titrino plus, 899 Coulometer.



1 pcs 6.2738.000 Funnel

For the KF coulometer.



1 pcs 6.2811.000 Molecular sieve

Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.



2 pcs 6.2816.030 Needle with Luer connector

Needle for filling in mercury into the Multi-Mode Electrode and for sample addition in Karl Fischer titration.



2 pcs **6.2816.090** **Syringe 5 mL**

5 mL syringe made of PP, with Luer connector.



Optional accessories

Order no.	Description
-----------	-------------

2.142.0100 Custom Q3X thermal printer

Compact printer with USB interface for

- 900 Touch Control
- 915 KF Ti-Touch
- 916 Ti-Touch
- 917 Coulometer
- 877 / 848 Titrino plus
- 865 / 876 Dosimat plus
- 91X Meter (cable 6.2151.140)
- Eco Dosimat / Titrator
- 862 Compact Titrosampler
- 870 KF Titrino plus
- 899 Coulometer



Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.

2.801.0010 801 Stirrer

Magnetic stirrer without stand for supplementing the Titrino plus, Dosimat plus, Titrandos, Sample Processors, 805 Dosimats and 780/781 pH Meters. With permanently attached cable for MSB (Metrohm Serial Bus).



2.803.0010 803 TI Stand with stirrer and pump

Compact titration stand with magnetic stirrer and built-in pump for solvent addition and aspiration of the titration vessel contents. Suitable for Karl Fischer titrations with Titrandos. MSB connection. Including comprehensive accessories, solvent bottle and waste bottle.



2.860.0010 860 KF Thermoprep

The 860 KF Thermoprep is designed for thermal sample preparation in Karl Fischer titration.

Numerous substances cannot be analysed by direct Karl Fischer titration as they are not soluble, react with the Karl Fischer reagent or release their water only very slowly or not until heated to high temperatures. Using tightly sealed sample vials, the sample is inserted in the oven. The samples can be subsequently analysed by any volumetric or coulometric KF titrator.



2.885.0010 885 Compact Oven Sample Changer

The 885 Compact Oven Sample Changer is used for thermal sample preparation in Karl Fischer titration. In the oven the samples can be heated up to 250°C. The moisture contained in the sample evaporates and is transported by a dry carrier gas into the titration cell, where the analysis takes place.



6.2141.390 Remote cable 885 - 899

Remote cable for connecting an 885 Compact Oven Sample Changer to the 899 Coulometer.



6.2148.030 RS-232/USB Box

Extends the Titrino plus by two RS-232 interfaces, e.g., for the connection of balances, PC (tiBase software), LIMS. Requires 6.2151.030 cable.



6.2164.500 Power Box for 899

The Power Box offers the option of operating the 899 Coulometer independently of the mains supply.



6.5406.000 Coulometric KF equipment, including generator electrode with diaphragm

Complete accessory set for KF coulometer.



6.5406.010 Coulometric KF equipment, including generator electrode without diaphragm

Complete accessory set for KF coulometer.

