


Accessories for A629281120

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery A629281120

Qty.	Order no.	Description
1 pcs	6.7440.480	NIRS Process, I/O Communication box
		I/O Communication box for process instruments with 24 analog and 24 digital channels



Optional accessories

Order no.	Description
A660044277	NIRS Analyzer PRO Lamp 12V/20W
Replacement lamp for the NIRS Analyzer Pro (12V/20W)	

