









Accessories for BWT-84000341

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery BWT-84000341

Qty.	Order no.	Description	
1 pcs	BWT-84000101	Laser safety goggles (785 nm)	
		Laser safety goggles for 785 nm lasers – suitable only for Class IIIB laser systems from B&W Tek.	
1 pcs	BWT-84000290	NanoRam Battery	
		NanoRam (model BWS456-785 only) Rechargeable lithium-ion battery.	
1 pcs	BWT-84000295	Vials for liquids, 15 mm in diameter	
		Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	BWT-84000296	Stylus/pen	
		Stylus/pen combo for touch screens.	
1 pcs	BWT-84000298	NanoRam Carrying Case	
		NanoRam hard plastic rugged carrying case.	
1 pcs	BWT-84000301	Power Adapter, 12 V DC	
		Power Adapter 12 V DC; 2.0 A.	
1 pcs	BWT-84000305	NanoRam ID Software	
		NID software for use on the PC for NanoRam account management, test result review, access to test reports and the Audit Trail.	

1 pcs	BWT-84000307	NanoRam Product Manual Set	
NanoRam documentation set including Product Manual and NanoRam ID Installation Instructions.			
1 pcs	BWT-84000329	Pouch for 15 mm Vials	
Pouch for 15 mm vials.			
1 pcs	BWT-84000343	Point-and-Shoot-Adapter	
Point-and-shoot adapter for handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.			
1 pcs	BWT-84000344	Vial Holder Adapter	
Vial holder adapter, compatible with 15 mm-diameter vials.			
1 pcs	BWT-84000345	NanoRam Bottle Adapter	
NanoRam Bottle Adapter.			
1 pcs	BWT-84000347	NanoRam Immersion Probe	
Immersion probe – 30.5 cm (12 inch) length – stainless steel with quartz objective.			
1 pcs	BWT-84000348	Polystyrene Validation Cap	
Polystyrene Validation cap with integrated ASTM standard polystyrene.			

1 pcs **BWT-84000352** **LAN cable**

LAN cable for connection to PC.



1 pcs **BWT-84000354** **Pouch for adapters**

Pouch for NanoRam adapters.



1 pcs **BWT-84000433** **O-ring replacement kit for NanoRam instruments**

NanoRam O-ring replacement kit – includes: O-rings, ID 9.00 mm, x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.



Optional accessories

Order no.	Description	
BWT-840000294	Battery Charging Cradle	
	Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
BWT-840000342	Right-Angle Adapter	
	Right-angle adapter for hands-free measurement using handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TactiClD-1064 models.	
BWT-840000346	NanoRam Tablet Holder	
	NanoRam tablet holder adapter.	
BWT-840000349	NanoRam 6 mm Shaft	
	NanoRam (model BWS456-785 only) replacement shaft with objective for sampling adapters – point-and-shoot; vial holder; bottle adapter. Working distance 6 mm.	
BWT-840000350	NanoRam 10 mm Shaft	
	NanoRam Shaft with lens for thick container sampling. Working distance 10 mm.	
BWT-840000353	Shaft Protection Cap	
	Handheld Raman shaft protective cap, package of 10.	
BWT-840000355	NanoRam Leather Shoulder Bag	
	Leather shoulder bag for NanoRam unit including one NR2-PCH-ADPT pouch for holding accessories.	

BWT-840000361 Handheld Laser Power Meter

Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.



BWT-840000481 NanoRam immersion shaft protective sleeve, package of 49

Disposable, protective sleeve for the NanoRam Immersion Probe model NR2-IMP, package of 49.



BWT-840000482 NanoRam immersion shaft protective sleeve, package of 490

Disposable, protective sleeve for the NanoRam Immersion Probe model NR2-IMP, package of 490.



BWT-840000655 Sterilized NanoRam Immersion Shaft Protective Sleeve

Sterilized disposable, protective sleeve for the NanoRam Immersion Probe model NR2-IMP; package of 490.



BWT-840000881 NanoRam Contact Probe

NanoRam accessories: 0.5" outer diameter x 150 mm L (5.9") – Raman contact probe accessory for NanoRam. Rugged 316 stainless steel shaft.

