

## Accessories for BWT-84000953

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

### Scope of delivery BWT-84000953

Qty.	Order no.	Description	
1 pcs	BWT-84000102	<b>Laser safety goggles (1064 nm)</b>  Laser safety goggles for lasers with excitation at 1,064 nm – suitable for Class IIIB laser systems from B&W Tek.	
1 pcs	BWT-84000294	<b>Battery Charging Cradle</b>  Lithium-ion battery charging cradle.(for charging NR-BAT and TID-BAT only).	
1 pcs	BWT-84000295	<b>Vials for liquids, 15 mm in diameter</b>  Borosilicate glass vials (diameter 15 mm), package of 6, with pouch.	
1 pcs	BWT-84000296	<b>Stylus/pen</b>  Stylus/pen combo for touch screens.	
1 pcs	BWT-84000298	<b>NanoRam Carrying Case</b>  NanoRam hard plastic rugged carrying case.	
1 pcs	BWT-84000329	<b>Pouch for 15 mm Vials</b>  Pouch for 15 mm vials.	
1 pcs	BWT-84000342	<b>Right-Angle Adapter</b>  Right-angle adapter for hands-free measurement using handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	

1 pcs	<b>BWT-84000343</b>	<b>Point-and-Shoot-Adapter</b>	Point-and-shoot adapter for handheld Raman instruments. Compatible with NanoRam, NanoRam-1064 and TacticID-1064 models.	
1 pcs	<b>BWT-84000344</b>	<b>Vial Holder Adapter</b>	Vial holder adapter, compatible with 15 mm-diameter vials.	
1 pcs	<b>BWT-84000348</b>	<b>Polystyrene Validation Cap</b>	Polystyrene Validation cap with integrated ASTM standard polystyrene.	
1 pcs	<b>BWT-84000354</b>	<b>Pouch for adapters</b>	Pouch for NanoRam adapters.	
1 pcs	<b>BWT-84000403</b>	<b>10.8 V rechargeable lithium ion battery</b>	Rechargeable lithium-ion battery, for use with NanoRam-1064, TacticID-GP Plus, TacticID-N Plus, TacticID-1064.	
1 pcs	<b>BWT-84000414</b>	<b>18 V Power Adapter</b>	Power adapter for 18 V DC, 1.67 A. Plugs for USA, EU, UK and AU included in the scope of delivery.	
1 pcs	<b>BWT-84000433</b>	<b>O-ring replacement kit for NanoRam instruments</b>	NanoRam O-ring replacement kit – includes: O-rings, ID 9.00 mm, x15; O-Rings, ID 10.00 mm x15; O-rings, ID 13.00 mm x10 for the shaft in handheld Raman instruments.	

---

**1 pcs    BWT-84000858    USB cable**

USB cable supporting micro USB port up to type A USB 2.0.



---

**1 pcs    BWT-84000969    Contact Immersion Probe (1064 nm)**

Handheld Raman Spectrometer 1064 accessories: Contact immersion probe, 6" length, working distance approx. 1 mm, operating pressure up to 2 bar at room temperature.



---

**1 pcs    BWT-84000970    NanoRam Large Bottle Adapter (1064)**

Handheld Raman Spectrometer 1064 accessories: Large bottle adapter with 10 mm working distance.




---

**1 pcs    BWT-84000971    NanoRam-1064 Software and Product Manual**

NanoRam-1064 accessories: USB flash drive containing NID EX software installation package; NanoRam-1064 Product Manual (PDF); NID EX Installation Instructions (PDF).



## Optional accessories

Order no.	Description	
<b>BWT-840000353</b>	<b>Shaft Protection Cap</b>	
	Handheld Raman shaft protective cap, package of 10.	
<b>BWT-840000361</b>	<b>Handheld Laser Power Meter</b>	
	Laser Power Meter, Pen Type, Calibrated. Range: 5 micro-Watts to 1 Watt. Used in OQ/PQ and annual certification laser power measurement.	