

## Accessories for 6.6072.209

Below, the accessories are grouped into Scope of delivery and Optional accessories.  
Please keep this printout at hand for ordering replacement material.  
These lists may be subject to change.

## Optional accessories

Order no.	Description
-----------	-------------

<b>6.6072.202</b>	<b>Vision Air 2.0 Pharma</b>
-------------------	------------------------------

### Vision Air - Universal spectroscopy software

Vision Air Pharma is a modern and simple-to-operate software solution for use in regulated environments.

Overview of the advantages of Vision Air Pharma:

- Individual software applications with adapted user interfaces ensure intuitive and simple operation
- Simple creation and maintenance of operating procedures
- SQL database for safe and simple data management
- In compliance with all requirements pursuant to Title 21 CFR Part 11

The Vision Air Pharma version (66072202) contains the following applications:

- Application for instrument and data management
- Application for routine analysis

Select this version when you already have a software for method development (Vision, PLS\_Toolbox or Unscrambler®) and work in regulated environments.

