





Accessories for 2.1001.1220

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1001.1220

Qty.	Order no.	Description	
1 pcs	1.1001.0020	OMNIS Titrator with stirrer	
		OMNIS Titrator with magnetic stirrer, without function license	
1 pcs	6.00404.300	dAg Titrode	
		Digital, combined silver ring electrode for OMNIS with a pH glass membrane as reference electrode. This maintenance-free electrode is suitable for precipitation titrations where the pH value remains constant (titrant silver nitrate), for example of <ul style="list-style-type: none"> • chloride, bromide, iodide • sulfides • hydrogen sulfide • mercaptans • cyanides 	
		This electrode is stored in distilled water.	
1 pcs	6.01600.010	OMNIS Liquid Adapter	
		Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.	
1 pcs	6.01601.000	Bottle cap, multi-use	
		Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.	

1 pcs	6.02004.010	Bottle holder	
			Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.
			
1 pcs	6.02005.000	OMNIS electrode holder	
			Electrode holder for one OMNIS titration system.
			
1 pcs	6.02005.010	Electrode Parking Station	
			Electrode Parking Station for electrode storage on an OMNIS titration system.
			
1 pcs	6.02100.010	Digital measuring module	
			Digital measurement channel for one OMNIS Titrator or Titration Module for the connection of digital "dTrodes" electrodes.
			
1 pcs	6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m	
			Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.
			
1 pcs	6.02107.000	Network cable, 1 m	
			Network cable (1 m) for connecting one OMNIS Sample Robot or OMNIS Titrator to an Ethernet Router, Ethernet Switch or local network
			

1 pcs 6.03001.220 OMNIS 20 mL cylinder unit

Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



1 pcs 6.06000.000 OMNIS Software storage device

License-based PC program for controlling complex OMNIS analysis systems.

- Sample-oriented operating procedures with numerous templates
- Graphical method editor with numerous templates
- Clearly organized and professional sample lists
- Management of instruments, work systems, sensors, solutions
- Dialog languages: German, English, Chinese



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.06001.910 Function license Advanced

"Advanced" function license for one OMNIS Titrator

- Liquid Handling (LQH)
- MEAS U/T/pH
- Calibration (Cal)
- Endpoint titration without conditioning (SET)
- Monotonic equivalence point titration (MET)
- Dynamic equivalence point titration (DET)



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.06003.010 OMNIS Stand-Alone license

Enables stand-alone operation of the OMNIS software on a Windows™ computer.

Features:

- The license already includes one OMNIS instrument license.
- Must be activated via the Metrohm licensing portal.
- Not transferable to another computer.



OMNIS
A WHOLE NEW LEVEL OF PERFORMANCE

1 pcs 6.1546.030 Extractor rod

For PTFE piston of dosing unit



1 pcs 6.1608.020 Glass bottle / 1000 mL / GL 45

Bottle for auxiliary solutions.



1 pcs 6.1608.023 Amber glass bottle / 1000 mL / GL 45

For Exchange Units. Bottle for auxiliary solutions.



1 pcs 6.1903.020 Stirring bar / 16 mm

Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 pcs 6.1903.030 Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



1 pcs 6.2013.010 Clamping ring

For support rods with a diameter of 10 mm.






1 pcs

6.2016.050

Support rod / 300 mm



Optional accessories

Order no.	Description	
6.05609.000	KFT equipment for OMNIS Titrator	
	Complete accessory set for volumetric Karl Fischer titration.	
6.05609.010	KFT equipment for OMNIS Titrator with Polytron	
	Complete accessory set for volumetric Karl Fischer titration with Polytron.	
6.06001.160	Function license KFT with conditioning	
	Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system.	