







Accessories for 2.1016.0110

Below, the accessories are grouped into Scope of delivery and Optional accessories.
Please keep this printout at hand for ordering replacement material.
These lists may be subject to change.

Scope of delivery 2.1016.0110

Qty.	Order no.	Description	
1 pcs	1.1016.0110	"Peristaltic" (4-channel) pump module	
		Module for installation in the module plate of the Pick&Place of the OMNIS Sample Robot. This workstation is equipped with two rinsing pumps and two aspiration pumps each. They are used to clean the sensors in two Pick&Place modules and to empty the sample beaker after the analysis before it is placed back in the rack once again.	
1 pcs	6.01803.000	Silicone tubing 8/12 2.5 m	
		Silicone tubing for connecting to the drainage of the Pick&Place module or pump module in order to channel any solution overflows that may occur safely into the waste canister.	
1 pcs	6.01804.500	Tubing adapter for drainage tubing	
		Tubing adapter for drainage tubing on Pick&Place modules and pump modules in the OMNIS Sample Robot Pick&Place.	
6 pcs	6.01805.000	FEP tubing connection / M6 / 66.5 cm	
		With light protection and protection against kinks.	
2 pcs	6.01805.010	PTFE tubing connection / M8 / 66.5 cm	
		With light protection and protection against kinks.	
2 pcs	6.05001.010	"Peristaltic" (2-channel) pump module accessory set	
		The accessory set contains all of the components required for the tube connections between a rinsing pump and an aspiration pump at the "Peristaltic" pump module. In addition to preassembled tubing with Luer connectors, each of these also includes two press clamps and pump tubing.	

2 pcs **6.1543.170** **Aspiration tip / M8 thread**

Aspiration tip for the Sample Processors.



4 pcs **6.1621.000** **Container 10 L**

As rinsing or waste container in automated systems.



4 pcs **6.1621.100** **Cap for 10-L container**

Cap for 10-L container. For easy and safe transportation of filled containers.



4 pcs **6.1828.000** **PVDF connection nipple**

For 6.1621.000 canister.



2 pcs **6.2709.070** **Guiding sleeve**

For fixation of buret and dosing tips in titration heads and electrode holders




6 pcs **6.2740.020** **Spray nozzle**

For the fine-spraying of the rinsing solution.



Optional accessories

Order no.	Description	
6.05700.010	Consumable Kit OMNIS Pump Module (2-channel)	
<p>Two pump tubings each and associated press clamps for replacement in the OMNIS pump module. One Consumable Kit can be used for re-equipping two peristaltic pumps of the pump module.</p>		
